











The Premium Sensor

Equipped with next-gen CMOS fiber optic technology, Clio sensors boast some of the highest quality images on the market, including the largest active area plates for maximum image coverage. The beveled edges and rounded corners of the Clio sensor's head give patients a more comfortable experience.

FEATURING

Easy Integration

Clio sensors seamlessly integrate with most dental imaging software. By integrating with your imaging software, our digital x-ray sensors are able to take advantage of your imaging software capabilities, making implementation and workflow easy.

Longer, Better Warranty & Replacement Program

2-Year Manufacturer Warranty

SOTA will repair or exchange parts or equipment due to manufacturer defect free of charge.

5-Year Complimentary Replacement Program Enrollment

Starting from the date of purchase, sensors that cannot be repaired and/or not covered under the 2-year manufacturer's warranty are offered significant discount pricing. Current rates are around 64% off MSRP price. (Contact SOTA for details).

Specs	Size 1	Size 2
External Dimensions	24.35 x 36.73 mm	30.49 x 42.80 mm
Active Area	19.95 x 30.02 mm	25.99 x 36.02 mm
Sensor Plate Thickness	5.41 mm	5.41 mm
Sensor Cable Length	2 meters	2 meters
Image Resolution	26.3 lp/mm theoretical	26.3 lp/mm theoretical
	20lp/mm actual	20lp/mm actual
Image Size	1580 x 1050 pixels	1896 x 1368 pixels
Image Pixel Size	19 x 19 microns	19 x 19 microns









* Available for other software platforms as well. Contact a SOTA representative for more details and availability.

Brand names, logos, and trademarks used herein remain the property of their respective owners. This listing of any firm or their logos is not intended to imply any endorsement or direct affiliation with SOTA Imaging.